

Human HSBP1 protein

Cat. No. GTX67479-pro

Species	Human
---------	-------

Package	100 µg
---------	--------

Applications

Observed MW (kDa)	9 kDa.
-------------------	--------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 50mM NaCl, 20% Glycerol
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	Full length protein, MAETDPKTVQ DLTSWQTL QQMQDKFQTM SDQIIGRIDD MSSRIDDLK NIADLMTQAG VEELESENKI PATQKS
-----------------	---

Expression System	E. coli
-------------------	---------

Purity	> 95% by SDS-PAGE.
--------	--------------------


Conjugation	Unconjugated
-------------	--------------

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX67479-pro Image

3µg Human HSBP1 protein (GTX67479-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 1