

Human HSF1 protein, His tag

Cat. No. GTX67481-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 57 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 50mM NaCl, 1mM DTT**Preservative** No preservative**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MDLPVGPAA GPSNVPAFLT KLWTLVSDPD
TDALICWSPS GNSFHVFDQG QFAKEVLPHY FKHNNMASFV RQLNMYGFRK WHIEQGGLV KPERDDTEFQ HPCFLRGQEQ
LLENIKRKVT SVSTLKSEDI KIRQDSVTKL LTDVQLMKGK QECMDSKLLA MKHENEALWR EVASLRQKHA QQQKVVNKLI
QFLISLVQSN RILGVKRIKIP LMLNDSGSAH SMPKYSRQFS LEHVHGSGPY SAPSPAYSSS SLYAPDAVAS SGPIISDITE
LAPASPMASP GGSIDERPLS SSPLVRVKEE PPSPQSPRV EEASPGRPSS VDTLLSPTAL IDSILRESEF APASVTALTD ARGHTDTEGR
PPSPPTSTP EKCLSVACLD KNELSDHLDA MDSNLDNLQT MLSSHGFSDV TSALLDLFSP SVTPDMSLP DLDSSLASIQ
ELLSPQEPPR PPEAENSSPD SGKQLVHYTA QPLFLDPGS VDTGSNDLPV LFELGEGSYF SEGDFGAEDP TISLLTGSEP PKAKDPTVS

Expression System E. coli**Purity** > 75% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67481-pro Image

3μg Human HSF1 protein (GTX67481-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).