

Human Hsp90 alpha protein, His tag

Cat. No. GTX67489-pro

Species Human**Package**

100 µg

Applications

Observed MW (kDa) 85 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 100mM NaCl**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MPEETQTQDQ PMEEEEVETF AFQAEIAQLM SLIINTFYSN KEIFLRELIS NSSDALDKIR YESLTDPSKL DSGKELHINL IPNKQDRTLIT IVDTGIGMTK ADLINNLGTI AKSGTKAFME ALQAGADISM IGQFGVGFYS AYLVAEKVTV ITKHNDDEQY AWESSAGGSF TVRTDTGEPM GRGTKVILHL KEDQTEYLEE RRIKEIVKKH SQFIGYPITL FVEKERDKEV SDDEAEEKED KEEKEKEEEK ESEDKPEIED VGSDEEEEEK DGDKKKKKKI KEKYIDQEEL NKTKEIPIWTRN PDDITNEEYG EFYKSLTNDW EDHLAVKHFS VEGQLEFRAL LFVPRRAPFD LFENRKKKNN IKLYVRRVFI MDNCEELIPE YLNFIRGVVD SEDLPLNISR EMLQQSKILK VIRKNLVKKC LELFTELAED KENYKKFYEQ FSKNIKLGIIH EDSQNRKKLS ELLRYYSAS GDEMVSCLKDY CTRMKENQKH IYYITGETKD QVANSAFVER LRKHGLEVIY MIEPIDEYCV QQLKEFEGKT LVSVTKEGLE LPEDEEEKKK QEEKKTKFEN LCKIMKDILE KKEKVVWSN RLVTSPPCIV TSTYGWTANM ERIMKAQALR DNSTMGYMAA KKHLEINPDH SIETLRQKA EADKNDKSVK DLVILLYETA LLSSGFSLEDPQTHANRIYR MIKLGGLIDE DDPTADDTSA AVTEEMPPLE GDDDDTSRMEE VD

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

DATA IMAGES



GTX67489-pro Image

3μg Human Hsp90 alpha protein (GTX67489-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).