

Human KPNA2 protein, His tag

Cat. No. GTX67510-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 58 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSHMSTNEN ANTPAARLHR FKNKGKDSTE MRRRRIEVNV ELRKAKKDDQ MLKRRNVSSF PDDATSPLQE NRNNQGTWNW SVDDIVKGIN SSNVENQLQA TQAARKLLSR EKQPPIDNII RAGLIPKFS FLGRTDCSPI QFESAWALTN IASGTSEQTK VVDGGAIPA FISLLASPHA HISEQAVWAL GNIAGDGSVF RDLVIKYGAV DPLLALLAVP DMSSLACGYL RNLTWTLNL CRNKNPAPPI DAVEQILPTL VRLHHDDPE VLADTCWAIS YLTGPNERI GMVKTGVWP QLVKLLGASE LPIVTPALRA IGNIVTGTDE QTQVVIDAGA LAVFPSLLTN PKTNIQKEAT WTMSNITAGR QDQIQQVNH GLVPFLSVL SKADFKTQKE AVWAVTNYTS GGTVEQIVYL VHCGIIEPLM NLLTAKDTKI ILVILDAISN IFQAAENLGE TEKLSIMIEE CGGLDKIEAL QNHENESVYK ASLSLIEKYF SVEEEDQNV VPETTSEGYT FQVQDGAPGT FNF

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67510-pro Image

3μg Human KPNA2 protein (GTX67510-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).