

Human SLP76 protein, His tag

Cat. No. GTX67519-pro

Species Human

Package
10 µg

Applications

Observed MW (kDa) 70 kDa.

Properties

Form Liquid

Buffer PBS, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, C-terminal His-Tag; ADPMALRNVP FRSEVLGWDP DSLADYFKKL NYKDCEKAVK KYHIDGARFL NLTENDIQKF PKLRVPILSK LSQINKNEE RRSIFTRKPQ VPRFPEETES HEEDNGGWSS FEEDDYESPND DDQDGEDDGD YESPNEEEEA PVEDDADYEP PPSNDEEALQ NSILPAKFPF NSNSMYIDRP PSGKTPQQPP VPPQRPMAAL PPPPAGRNHS PLPPPQTNHE EPSRSRNHKT AKLPAPSIDR STKPPLDRSL APFDREPFTL GKKPPFSDKP SIPAGRSLGE HLPKIQKPPL PPTTERHERS SPLPGKKPPV PKHGWGPD RR ENEDDDVHQR PLQPALLPM SSNTFPSRST KPSPMNPLPS SHMPGAFSES NSSFPQSASL PPYFSQGGSN RPIRAEGRN FPLPLPNKPR PPSPAE EENS LNEEWYVSYI TRPEAEAALR KINQDGTFLV RDSSKKTITN PYVLMVLYKD KVYNIQIRYQ KESQVYLLGT GLRGKEDFLS VSDIIDYFRK MPLLIDGKN RGSRYQCTLT HAAGYPHHHH HH

Expression System Baculovirus

Purity > 90% by SDS-PAGE.

Endotoxin < 1.0 EU/µg (determined by LAL method)

Conjugation Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES

GTX67519-pro Image

3μg Human SLP76 protein (GTX67519-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).