

Human LTA4H protein, His tag

Cat. No. GTX67525-pro

Species Human**Package**

20 µg

Applications

Observed MW (kDa) 69 kDa.

Properties

Form Liquid**Buffer** PBS, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPDIVDT CSLASPASVC RTKHLHLRCS
VDFTRRTLGT TAALTQSQE DNLRLSLVDT KDLTIKWI NGQEVKYGALG ERQSYKGSPM EISLPALSK NQEIIEISF ETSPKSSALQ
WLTPEQTSKG EHPYLFSCQ AIHCRAILPC QDTPSVKLT TAEVSVPKEL VALMSAIRDG ETPDPEDPSR KIYKFIQKVP IPCYLIALV
GALESRQIGP RTLWSEKEQ VEKSAIEFSE TESMLKIAED LGGPYVWGQY DLLVLPSPF YGGMENPCLT FVTPTLLAGD
KSLSNVIAHE ISHSWTGNLV TNKTWDHFWL NEGHTVYLER HICGRLEGEK FRHFNALGGW GELQNSVKTF GETHPFTKLV
VDLTDIDPDV AYSSVPYEKG FALLFYLEQL LGGPEIFLGF LKAYVEKFSY KSITDDWWD FLYSYFKDKV DVLNQVDWNA
WLYSPGLPPI KPNYDMLTN ACIALSQRWI TAKEDDLNSF NATDLKDLSS HQLNEFLAQT LQRAPLPLGH IKRMQEVYNF
NAINNSEIRF RWLRLCIQSK WEDAIPLALK MATEQGRMKF TRPLFKDLAA FDKSHDQAVR TYQEHKASMH PVTAMLVGKD LKVD

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67525-pro Image

3μg Human LTA4H protein (GTX67525-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).