

Human Myosin Light Chain 2 (MLC-2v) protein, His tag

Cat. No. GTX67569-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 19 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 40% Glycerol, 5mM CaCl₂.

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MAPKKAKKRA GGANSNVFSM FEQTQIQEFK
EAFTIMDQNR DGFIDKNDLR DTFAALGRVN VKNEEIDEMI KEAPGPINFT VFLTMFGEKL KGADPEETIL NAFKVFDPPEG
KGVLKADYVR EMLTTQAERF SKEEVDQMFA AFPPDVTGNL DYKNLVHIIT HGEEKD

Expression System E. coli

Purity > 95% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67569-pro Image

3μg Human Myosin Light Chain 2 (MLC-2v) protein (GTX67569-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).