

Human NDUFA2 protein, His tag

Cat. No. GTX67580-pro

Species Human**Package**

50 µg

Applications


Observed MW (kDa) 11 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAAAAS RGVGAKLGLR EIRIHLQRS
PGSQGVRDFI EKRYVELKKA NPDLPIIRE CSDVQPKLWA RYAFGQETNV PLNNFSADQV TRALENVLSG KA**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX67580-pro Image

3µg Human NDUFA2 protein (GTX67580-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 08 Page 1 of 1