

# Human NMT1 protein, His tag

**Cat. No. GTX67597-pro**
**Species** Human

**Package**  
100 µg

**Applications**
**Observed MW (kDa)** 57 kDa.

**Properties**
**Form** Liquid

**Buffer** PBS

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMADSEET AVKPPAPPLP QMMEGNGNGH EHCSDCENEE DNSYNRGGLS PANDTGAKKK KKKQKKKKEK GSETDSAQDQ PVKMNSLPAE RIQEIQKAIE LFSVGQGPAK TMEEASKRSY QFWDTPVPK LGEVNTGHP VEPDKDNIRQ EPYTLPGGFT WDALDLGDRG VLKELYTLN ENYVEDDDNM FRFDYSPEFL LWALRPPGWL PQWHCGVRV SSRKLVGFIS AIPANIHIYD TEKKMVEINF LCVHKKLRSK RVAPVLIREI TRRVHLEGIF QAVYTAGVWL PKPVGTCRYW HRSLNPRKLI EVKFSHLRN MTMQRTMKLY RLPETPKTAG LRPMETKDIP VHQLLTRYL KQFHLTPVMS QEEVEHWFYP QENIIDTFV ENANGEVTDF LSFYTLPSI MNHPTHKSLK AAYSFYNVHT QTPLLDLMSD ALVLAKMKGF DVFNALDLME NKTFLKLFK GIGDGNLQYY LYNWKCPMSG AEKVGLVLQ

**Expression System** E. coli

**Purity** > 95% by SDS-PAGE.

**Conjugation** Unconjugated

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX67597-pro Image

3μg Human NMT1 protein (GTX67597-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).