

Human PDE6H protein, His tag

Cat. No. GTX67621-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 9 kDa.


Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 40% Glycerol, 2mM DTT, 0.1mM PMSF**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGMSDNTTL PAPASNQGPT TPRKGPPKFK QRQTRQFKSK PPKKGVKGF G DDIPGMEGLG TDITVICPWE AFSHLELHEL AQFGII**Expression System** E. coli**Purity** > 80% by SDS-PAGE.**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX67621-pro Image**

3µg Human PDE6H protein (GTX67621-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

For full product information, images and publications, please visit our [website](https://www.genetex.com).