

Human PEPD protein, His tag

Cat. No. GTX67623-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 55 kDa.

Properties

Form Liquid

Buffer PBS, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAATGP SFWLGNETLK VPLALFALNR QRLCERLRKN PAVQAGSIV LQGGEETQRY CTDGVLFRQ ESFFHWAQGV TEPGCYGVID VDTGKSTLFV PRLPASHATW MGKIHSEHF KEKYAVDDVQ YVDEIASVLT SQKPSVLLTL RGVNTDSGSV CREASFDGIS KFEVNTILH PEIVECRVFK TDMELEVLRY TNKISSEAHR EVMKAVKVGM KEYEESLFE HCYSRGGMR HSSYTCICGS GENSAVLHYG HAGAPNDRTI QNGDMCLFDM GGEYYCFASD ITCSPFANGK FTADQKAVYE AVLRSSRAVM GAMKPGVWWP DMHRLADRIH LEELAHMGIL SGSVDAMVQA HLGAVFMPHG LGHFLGIDVH DVGGYPEGVE RIDEPLRLSL RTARHLQPGM VLTVEPGIYF IDHLLDEALA DPARASFLNR EVLQRFGRFG GVRIEEDVVV TDSGIELLTC VPRTVEEIEA CMAGCDKAPT PFSGPK

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67623-pro Image

3μg Human PEPD protein (GTX67623-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).