

Human PGM1 protein, His tag

Cat. No. GTX67635-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 61 kDa.

Properties

Form Liquid

Buffer PBS, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMVKIVTV KTQAYQDQKP GTSGLRKRVK VEQSSANYAE NFIQSIISTV EPAQRQEATL VGGDGRFYM KEAQLIARI AAANGIGRLV IGQNGILSTP AVSCIIRKIK AIGGIILTAS HNPGGPNGDF GIKFNISNGG PAPEAITDKI FQISKTEEY AVCPDLKVDL GVLGKQQFDL ENKFKPFTVE IVDSVEAYAT MLRSIFDFSA LKELLSGPNR LKIRIDAMHG VGPYVKKIL CEELGAPANS AVNCVPLEDF GGHHPDPNLT YAADLVETMK SGEHDFGAFF DGDGDRNMIL GKHGFFVNPS DSVAVIAANI FSIPLYQQTG VRGFARSMPT SGALDRVASA TKIALYETPT GWKFFGNLMD ASKLSLCGEE SFGTGSDBIR EKDGLWAVLA WLSILATRKQ SVEDILKDHQ QKYGRNFFTR YDYEEVEAEG ANKMMKDLEA LMFDRSFVGK QFSANDKVYT VEKADNFEYS DPVDGSSIRN QGLRLIFTDG SRIVFRLSGT GSAGATIRLY IDSYEKDVAK INQDPQVMLA PLISIALKVS QLQERTGRTA PTVIT

Expression System E. coli

Purity > 95% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67635-pro Image

3μg Human PGM1 protein (GTX67635-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).