

Human POLR2K protein, His tag

Cat. No. GTX67654-pro
Species Human

Package
50 µg

Applications
Observed MW (kDa) 7 kDa.

Properties
Form Liquid

Buffer PBS, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMDTQKDV QPPKQQPMIY ICGECHTENE IKSRDPIRCR ECGYRIMYKK RTKRLVFDA R

Expression System E. coli


Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX67654-pro Image

3µg Human POLR2K protein (GTX67654-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).