

Human POU6F1 protein, His tag

Cat. No. GTX67657-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 33 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPGISSQ ILTNAQGQVI GTLPWWVNSA
SVAAPAPAQS LQVQAVTPQL LLNAQGQVIA TLASSPLPPP VAVRKPSTPE SPAKSEVQPI QPTPTVPQPA VVIASPAPAA KPSASAPIPI
TCSETPTVSQ LVSKPHTPSL DEDGINLEEI REFAKNFKIR RLSLGLTQTQ VGQALTATEG PAYQSAICR FEKLDITPKS AQKLKPVLEK
WLNEAELRNQ EGQQNLMEFV GGEPSSKKRRK RTSFTPQAIE ALNAYFEKNP LPTGQEITEI AKELNYDREV VRWWFCNRRQ
TLKNTSKLNV FQIP

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67657-pro Image

3µg Human POU6F1 protein (GTX67657-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).