

Human Calcineurin A protein, His tag

Cat. No. GTX67664-pro

Species Human**Package**

10 µg

Applications

Observed MW (kDa) 59 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 20% Glycerol, 5mM DTT, 1mM EDTA**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH TGSMSSEPKAI DPKLSTTDRV VKAVPFPPSH
RLTAKEVFDN DGKPRVDILK AHLMKEGRL ESVALRIITE GASILRQEK LLDIDAPVTV CGDIHGQFFD LMKLFEVGGG PANTRYLFLG
DYVDRGYFSI ECVLYLWALK ILYPKTLFLL RGNHECRHLT EYFTFKQECK IKYSERVYDA CMDAFDCLPL AALMNQQFLC
VHGGLSPEIN TLDDIRKLDR FKEPPAYGPM CDILWSDPLE DFGNEKTQEH FTHNTVRGCS YFYSYPAVCE FLQHNNLLSI
LRAHEAQDAG YRMYRKSQTT GFPSLITFS APNYLDVYNN KAAVLKYENN VMNIRQFNCS PHPYWLPNFM DVFTWSLPFV
GEKVTEMLVN VLNICSDDEL GSEEDGFDGA TAAARKEVIR NKIRAIGKMA RVFSLREES ESVLTLKGLT PTGMLPSGVL
SGGKQTLQSA IKGFSPOHKI TSFEEAKGLD RINERMPPRR DAMPSDANLN SINKALTSET NGTDSNGSNS SNIQ

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67664-pro Image

3μg Human Calcineurin A protein (GTX67664-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).