

Human JNK3 protein, His tag

Cat. No. GTX67680-pro

Species Human**Package**

10 µg

Applications

Observed MW (kDa) 53 kDa.

Properties

Form Liquid**Buffer** PBS, 30% Glycerol, 1mM DTT, 0.1mM PMSF**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, C-terminal His-Tag; MSLHFLYYCS EPTLDVKIAF CQGFDKQVDV SYIAKHYNMS KSKVDNQFYS
VEVG DSTFTV LKRYQNLKPI GSGAQGIVCA AYDAVLDRNV AIKKLSRPFQ NQTHAKRAYR ELVLMKCVNH KNIISLLNVF
TPQKTL EEFQ DVYLVME LMD ANLCQVIQME LDHERMSYLL YQMLCGIKHL HSAGIIHRDL KPSNIVVKSD CTLKILDFGL
ARTAGTSFMM TPYVVTRYR APEVILMGY KENVDIWSVG CIMGEMVRHK ILFPGRDYID QWNKVIEQLG TPCPEFMKKL
QPTVRNYVEN RPKYAGLTFF KLPDLSLFA DSEHNKLKAS QARDLLSKML VIDPAKRISV DDALQHPYIN WYDPAEVEA
PPPQIYDKQL DEREHTIEEW KELIYKEVMN SEEKTKNGVV KGQPSPSGAA VNSSESLPPS SSVNDISSMS TDQTLASDTD
SSLEASAGPL GCCRHHHHHHH

Expression System Baculovirus**Purity** > 85% by SDS-PAGE.**Endotoxin** < 0.1 EU/µg (determined by LAL method)**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67680-pro Image

3μg Human JNK3 protein (GTX67680-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).