

Human Kallikrein 7 protein, His tag

Cat. No. GTX67688-pro

Species	Human
---------	-------

Package	20 µg
---------	-------

Applications

Observed MW (kDa)	30 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	PBS, 10% Glycerol
--------	-------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	----------------------------------------------------------------------------

Region/Sequence	Full length protein, C-terminal His-Tag; ADPMNEYTVH LGSDTLGDRR AQRIKASKSF RHPGYSTQTH VNDLMLVKLN SQARLSSMVK KVRLPSCRCEP PGTCTVSGW GTTSPDVTF PSDLMCV DVK LISPQDCTKV YKDLENSML CAGIPDSKKN ACNGDSGGPL VCRGTLQGLV SWGTFCGQP NDPGVYTQVC KFTKWINDTM KKHRHHHHHH
-----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Expression System	Baculovirus
-------------------	-------------

Purity	> 95% by SDS-PAGE.
--------	--------------------

Endotoxin	< 1.0 EU/µg (determined by LAL method)
-----------	----------------------------------------

Conjugation	Unconjugated
-------------	--------------


Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---------------------------------------------------------------------------------------------------------------------------------------------------

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67688-pro Image

3μg Human Kallikrein 7 protein (GTX67688-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).