

Human SHP2 protein, His tag

Cat. No. GTX67711-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 54 kDa.

Properties

Form Liquid**Buffer** PBS, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, C-terminal His-Tag; ADPMTSRWF HPNITGVEAE NLLLTRGVDG SFLARPSKSN PGDFTLSVRR
NGAVTHIKIQ NTGDYYDLYG GEKFATLAEL VQYMEHHGQ LKEKNGDVIE LKYPLNCADP TSEWTFHGH LSGKEAEKLLT
EKGKHGSFLV RESQSHPGDF VLSVRTGDDK GESNDGKSKV THVMIRCQEL KYDVGGGERF DSLTDLVEHY KKNPMVETLG
TVLQKQPLN TTRINAAEIE SRVRELSKLA ETDDKVKQGF WEEFETLQQQ ECKLLYSRKE GQRQENKNKN RYKNILPFDH
TRVVLHDGDP NEPVSDYINA NIIMPEFETK CNNSKPKKSY IATQGCLQNT VNDFWRMV FQ ENSRVIVMTT KEVERGKSKC
VKYWPDEYAL KEYGVMRVRN VKESAHDYT LRELKLSKVG QGNTERTVWQ YHFRTWPDHG VPSPGGVLD FLEEVHHKQE
SIMDAGPVVV HCRHHHHHHH

Expression System Baculovirus**Purity** > 90% by SDS-PAGE.**Endotoxin** < 1.0 EU/µg (determined by LAL method)**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67711-pro Image

3μg Human SHP2 protein (GTX67711-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).