

Human PYGL protein, His tag

Cat. No. GTX67717-pro

Species Human**Package**

10 µg

Applications

Observed MW (kDa) 97 kDa.

Properties

Form Liquid**Buffer** PBS, 30% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSEFELRRQ ASMAKPLTDQ EKRRQISIRG
IVGVENVAEL KKSFNRLHF TLVKDRNVAT TRDYFALAH TVRDHLVGRW IRTQQHYDYK CPKRYYLSL EFYMGRTLQN
TMINLGLQNA CDEAIYQLGL DIEELEIEE DAGLGNGGLG RLAACFLDSM ATGLAAYGY GIRYEGIFN QKIRDGWQVE
EADDWLRYN PWEKSRPEFM LPVHFYGVKE HTNTGTKWID TQWLALPYD TPVPGYMNNT VNTMRLWSAR APNDFNLRDF
NVGDYIQAVL DRNLAENISR VLYPNDNFFE GKELRLKQEY FVVAATLQDI IRRFKASKFG STRGAGTVFD AFPDQVAIQL
NDTHPALAIP ELMRIFVDIE KLPWSKAWEL TQKTFAYTNH TVLPEALERW PVDLVEKLLP RHLEIYEIN QKHLDRIVAL FPKDVDRLRR
MSLIEEGSK RINMAHLCIV GSHAVNGVAK IHSDIVKTKV FKDFSELEPD KFQNKTNGIT PRRWLLCNP GLAELIAEKI GEDYVKDLSQ
LTKLSFLGD DVFLRELAKV KQENKLKFSQ FLETEYKVKI NPSSMFDVQV KRIHEYKRQL LNCLHVITMY NRIKKDPKKL FVPRTVIIGG
KAAPGYHMAK MIIKLITSA DVVNNDPMVG SKLKVIFLEN YRVSLAEKVI PATDLSEQIS TAGTEASGTG NMKFMLNGAL
TIGTMDGANV EMAGEEAGEEN LFIFGMRIDD VAALDKKGYE AKEYYEALPE LKLVIDQIDN GFFSPKQPD LFKDIINMLFY
HDRFKVFADY EAYVKCQDKV SQLYMNP KAW NTMVLKNIAA SGKFSSDRTI KEYAQNIWNV EPSDLKISLS NESNKVNGN

Expression System E. coli**Purity** > 80% by SDS-PAGE.**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES



GTX67717-pro Image

3μg Human PYGL protein (GTX67717-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).