

Human QDPR protein, His tag

Cat. No. GTX67718-pro

Species Human

Package

100 µg

Applications

Observed MW (kDa) 26 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAAAAA GEARRLVYVY GRGALGSRCV
 QAFRRARNWWW ASVDWVENE ASASIIVKMT DSFTEQADQV TAEVGKLLGE EKVDAILCVA GGWAGGNAKS KSLFKNCDLM
 WKQSIWTSTI SSSLATKHLK EGGLTLAGA KAALDGTGPM IGYGMAKGAV HQLCQSLAGK NSGMPPGAAA IAVLPVTLDT
 PMNRKSMPEA DFSSWTPLEF LVETFHDWIT GKNRPSSGSL IQVWTEGRT ELTPAYF

Expression System E. coli

Purity > 90% by SDS-PAGE.


Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX67718-pro Image

3µg Human QDPR protein (GTX67718-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).