

Human CRALBP protein, His tag

Cat. No. GTX67752-pro
Species Human

Package
50 µg

Applications
Observed MW (kDa) 40 kDa.

Properties
Form Liquid

Buffer PBS, 20% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, C-terminal His-Tag; ADPMSEGVGT FRMPVEEQE LRAQLEQLTT KDHGPVFGPC SQLPRHTLQK AKDELNEREE TREEAVRELQ EMVQAQAASG EELAVAVAER VQEKDSGFFL RFIRARKFNV GRAYELLRGY VNFRLQYPEL FDSLSPFAVR CTIEAGYPGV LSSRDYGRV VMLFNINWQ SQEITFDEIL QAYCFLEKL LENEETQING FCIENFKGF TMQQAASLRT SDRKMMVDML QDSFPAFKA IHFIHQPWYF TTYNVVKPF LKSKLLERVF VHGDLSGFY QEIDENILPS DFGGTLPKYD GKAVAEQLFG PQAQAFNTAF HHHHHH

Expression System Baculovirus

Purity > 90% by SDS-PAGE.

Endotoxin < 1.0 EU/µg (determined by LAL method)

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67752-pro Image

3μg Human CRALBP protein (GTX67752-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).