

# Human RPL18A protein, His tag

**Cat. No. GTX67761-pro**

**Species** Human

**Package**  
100 µg

## Applications

**Observed MW (kDa)** 21 kDa.

## Properties

**Form** Liquid

**Buffer** 20mM Tris-HCl, 0.4M Urea, 10% Glycerol

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence** Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMKASGTL REYKVWGRCL PTPKCHTPPL YRMRIAPNH VAKSRFWYF VSQKKMKKS SGEIVYCGQV FEKSPLRVKN FGIWLYDSR SGTHNMYREY RDLTTAGAVT QCYRDMGARH RARAHSIQIM KVEEIAASKC RRAVKQFHD SKIKFPLPHR VLRRQHKPRF TTKRPNTFF

**Expression System** E. coli

**Purity** > 85% by SDS-PAGE.


**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX67761-pro Image**

3μg Human RPL18A protein (GTX67761-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).