

Human RPS14 protein, His tag

Cat. No. GTX67780-pro

Species Human

Package

10 µg

Applications

Observed MW (kDa) 16 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 40% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAPRK GK EKKEEQVISL GPQVAEGENV FGVCHIFASF NDTFVHVTDL SGKETICRVT GGMKVKADRD ESSPYAAMLA AQDVAQRCKE LGITALHIKL RATGGNRTKT PGPGAQSALR ALARSGMKIG RIEDVTPIPS DSTRRKGGRR GRRL

Expression System E. coli

Purity > 80% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67780-pro Image

3μg Human RPS14 protein (GTX67780-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).