

Human S100A3 protein, His tag

Cat. No. GTX67788-pro

Species Human**Package**
100 µg

Applications


Observed MW (kDa) 12 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 20% Glycerol, 2mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MARPLEQAVA AIVCTFQEYA GRCGDKYKLC QAEKELLQK ELATWTPTEF RECDYNKFMS VLDTNKDCEV DFVEYVRS LA CLCLYCHEYF KDCPSEPPCS Q**Expression System** E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX67788-pro Image

3µg Human S100A3 protein (GTX67788-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 07 Page 1 of 1