

Human alpha Synuclein protein

Cat. No. GTX67834-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 14 kDa.


Properties

Form Lyophilized powder**Buffer** Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 0.2µm-filtered 20mM Tris-HCl buffer (pH7.5), 100mM NaCl, 1mM MgCl₂.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** Full length protein, MDVFMKGLSK AKEGVAAAE KTKQGVAAEA GKTKEGVLYV GSKTKEGVVH GVATVAEKT EQVTNVGGAV VTGVTAVAQK TVEGAGSIAA ATGFVKDQL GKNEEGAPQE GILEDMPVDP DNEAYEMPSE EGYQDYEPEA**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Endotoxin** < 1.0 EU/µg (determined by LAL method)**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX67834-pro Image**

3µg Human alpha Synuclein protein (GTX67834-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

For full product information, images and publications, please visit our [website](https://www.genetex.com).