

Human SNRPA1 protein, His tag

Cat. No. GTX67837-pro

Species Human

Package
10 µg

Applications

Observed MW (kDa) 28 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 200mM NaCl, 30% Glycerol, 2mM DTT, 0.1mM PMSF

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MVKLTAEIE QAAQYTNVR DRELDLRGYK IPVIENLGAT LDQFDAIDFS DNEIRKLDGF PLLRRLKTL VNNNRICRIG EGLDQALPCL TELILTNNLSL VELGDLDPLA SLKSLTYLSI LRNPVTNKKH YRLYVIYKVP QVRVLDQKLV KLKERQAEK MFKGKRGAL AKDIARRSKT FNPAGALPTD KKKGGPSPGD VEAIKNAIAN ASTLAEVERL KGLLQSGQIP GRERRSGPTD DGEEMEEDT VTNGS

Expression System E. coli

Purity > 80% by SDS-PAGE.

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67837-pro Image

3μg Human SNRPA1 protein (GTX67837-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).