

Human Aurora A protein, His tag

Cat. No. GTX67859-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 46 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, 0.5mM DTT, 0.1mM EDTA, 0.1mM EGTA, 0.1mM PMSF**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MDRSKENCIS GPVKATAPVG GPKRVLVTQQ
FPCQNPLPVN SGQAQRVLCF SNSSQRIPLQ AQKLVSHPK VQNQKQKQLQ ATSVPHVSR PLNNTQKSKQ PLPSAPENNP
EEELASKQKN EESKKRQWAL EDFEIGRPLG KGKFGNVYLA REKQSKFILA LKVLFAQLE KAGVEHQLRR EVEIQSHLRH PNILRLYGIF
HDATRVYLIL EYAPLGTVYR ELQKLSKFDE QRTATYITEL ANALSYCHSK RVIHRDIKPE NLLLSAGEL KIADFGWSVH APSSRRITLC
GTLDYLPPEM IEGRMHDEKV DLWSLGVLCY EFLVGKPPFE ANTYQETYKR ISRVEFTFPD FVTEGARDLI SRLLKHNPQS
RPMLREVLEH PWITANSSKP SNCQNKESAS KQS

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67859-pro Image

3μg Human Aurora A protein (GTX67859-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).