

Human TRBP2 protein

Cat. No. GTX67869-pro

Species	Human
---------	-------

Package	100 µg
---------	--------

Applications

Observed MW (kDa)	39 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 0.4M Urea, 10% Glycerol
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	Full length protein, MLAANPGKTP ISLLQEYGTR IGKTPVYDLL KAEGQAHQPN FTFRVTVGDT SCTGQGPSKK AAKHKAAEVA LKHLKGGSML EPALEDSSSF SPLDSSLPED IPVFTAAAAA TPVPSVVLTR SPPMELQPPV SPQQSECNPV GALQELVWQK GWRLPEYTVT QESGPAHRKE FTMTCTRVERF IEIGSGTSKK LAKRNAAAKM LLRVHTVPLD ARDGNVEPD DDHFSIGVGS RLDGLRNRP GCTWDSLRNS VGEKILSLRS CSLGSLGALG PACCRVSEL SEEQAFHVS YLDIEELSLSG LCQCLVELST QPATVCHGSA TTREAARGEA ARRALQYLKI MAGSK
-----------------	--

Expression System	E. coli
-------------------	---------

Purity	> 85% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------


Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67869-pro Image

3μg Human TRBP2 protein (GTX67869-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).