

Human ThrRS protein, His tag

Cat. No. GTX67870-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 83 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 2mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MFEEKASSPS GKMGGEEKPI GAGEEKQKEG GKKNKEGSG DGGRAELNPW PEIYTRLEM YNILKAHDS ILAEKAEKDS KPIKVTLPDG KQVDAESWKT TPYQIACGIS QGLADNTVIA KVVNNVWDL D RPLEEDCTLE LLKFEDEEAQ AVYWHSSAHI MGEAMERVY GCLCYGPPIE NGFYDMMYLE EGGVSSNDFS SLEALCKKII KEKQAFERLE VKKETLLAMF KYNKFKCRIL NEKVNTPTTT VYRCGPLIDL CRGPHVRHTG KIKALKIHKN SSTYWEKGAD METLQRIYGI SFPDPKMLKE WEKFQEEAKN RDHRKIGRDQ ELYFFHELSP GSCFFLPKGA YIYNALIEFI RSEYRKRGFQ EVVTPNIFNS RLWMTSGHWQ HYSNMFSFE VEKELFALKP MNCPGHCLMF DHRPSWREL PLRLADFGVL HRNELSGALT GLTRVRRFQQ DDAHIFCAME QIEDEIKGCL DFLRTVYSVF GFSFKLNLST RPEKFLGDIE VWDQAEKQLE NSLNEFGEKW ELNSGDGAFY GPKIDIQKD AIGRYHQCAT IQLDFQLPIR FNLTVVSHDG DDKKRPVIVH RAILGSVERM IAILTENYGG KWPFWLSPRQ VMVVPVGPTC DEYAQKVRQQ FHDAKFMADI DLDPGCTLNK KIRNAQLAQY NFILWGEKE KISGTVNI RT RDNKVHGERT ISETIERLQQ LKEFRSKQAE EEF

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 09 / 07 Page 1 of 2

DATA IMAGES



GTX67870-pro Image

3μg Human ThrRS protein (GTX67870-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).