

Human TDG protein, His tag

Cat. No. GTX67882-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 46 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 400mM Urea, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMEAENAG SYSLQQAQAF YTFPFQQLMA
EAPNMAVVNE QQMPPEVPAP APAQEPVQEA PKGRKRKPRT TEPKQPVEPK KPVESKKS GK SAKSKEKQEK ITDTFKVKRK
VDRFNGVSEA ELLTKTLPDI LTFNLDIVII GINPGLMAAY KGHYPGPGN HFWKCLFMSG LSEVQLNHMD DHTLPGKYGI
GFTNMVERTT PGSKDLSSKE FREGGRILVQ KLQKYQPRIA VFNGKCIYEI FSKEVFGVKV KNLEFGLQPH KIPDTETLCY
GMPSSSARCA QFPRAQDKVH YYIKLDLRD QLKGIERNMD VQEVQYTFDL QLAQEDAKKM AVKEEKYDPG YEAAYGGAYG
ENPCSSEPCG FSSNGLIESV ELRGESAFSG IPNGQWMTQS FTDQIPFSFN HCGTQEQQEE SHA

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx67882-pro Image

3μg Human TDG protein (GTx67882-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).