

Human TUBG1 protein, His tag

Cat. No. GTX67910-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 52 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 100mM NaCl, 40% Glycerol, 2mM DTT, 50mM Imidazole**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, C-terminal His-Tag; MPREITLQL GQCGNQIGFE FWKQLCAEHG ISPEGIVEEF ATEGTRKDV FFYQADDEHY
IPRAVLDDLE PRVIHSILNS PYAKLYNPEN IYLSEHGGA GNNWASGFSQ GEKIHEDIFD IIDREADGSD SLEGFVLCHS IAGGTGSGLG
SYLLERLNDR YPKKLVTYS VFPNQDEMSD VVQPYNSLL TLKRLTNAD CVWLDNTAL NRIATDRLHI QNPSFSQINQ
LVSTIMSAST TTLRYPGYMN NDLIGLIASL IPTPRHLFLM TGYTPLTDDQ SVASVRKTTV LDVMRRLLQP KNVMVSTGRD
RQTNHCYIAI LNIIQGEVDP TQVHKSLQRI RERKLANFIP WGPASIQVAL SRKSPYLP SA HRVSGLM MAN HTSISSLFER
TCRQYDKLRK REAFLEQFRK EDMFKDNFDE MDTSREIVQQ LIDEYHAATR PDYISWGTQE QHHHHHHH

Expression System Baculovirus**Purity** > 85% by SDS-PAGE.**Endotoxin** < 1.0 EU/µg (determined by LAL method)**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67910-pro Image

3μg Human TUBG1 protein (GTX67910-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).