

Human UGP2 protein, His tag

Cat. No. GTX67928-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 57 kDa.

Properties

Form Liquid**Buffer** PBS, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMSRFVQD LSKAMSQDGA SQFQEVIRQE
LELSVKKELE KILTASSHE FEHTKKDLTG FRKLFHRFLQ EKGPSVDWVGK IQRPPEDSIQ PYEKIKARGL PDNISSVLNK LVVVKLNGGL
GTSMGCKGPK SLIGVRNENT FLDLTVQQIE HLNKTYNTDV PLVLMNSFNT DEDTKKILQK YNHCRVKIYT FNQSRYPRI
KESLLPAKDV SYSGENTEA WYPPGHGDIY ASFYNSGLLD TFIGEGKEYI FVSNIDNLGA TVDLYLNLHL MNPPNGKRCE
FVMEVTNKTR ADVKGGTLTQ YEGKLRLVEI AQVPAHVDE FKSVSFKIF NTNNLWISLA AVKRLQEQA IDMEIIVNAK
TLDGGLNVIQ LETAVGAAIK SFENSLGINV PRSRFLPVKT TSDLLVMSN LYSLNAGSLT MSEKREFPTV PLVKLGSSFT KVQDYLRRE
SIPDMEILDH LTVSGDVTFG KNVSLKGTVI IIANHGDRID IPPGAVLENK IVSGNLRILD H

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67928-pro Image

3μg Human UGP2 protein (GTX67928-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).