

Human ZFAND5 protein, His tag

Cat. No. GTX67950-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 23 kDa.

Properties

Form Liquid

Buffer PBS, 30% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMQETNQ TPGPMLCSTG CGFYGNPRTN GMCSVCKYKEH LQRQQNSGRM SPMGTASGSN SPTSDSASVQ RADTSLNNCE GAAGSTSEKS RNPVPAALPV TQQMTEMESIS REDKITTPKT EVSEPVWTQP SPSVSQPSTS QSEKAPLKP KKKKNRCFMC RKKVGLTGFD CRCGNLFCGL HRYSDKHNCY YDYKAEAAAK IRKENPVVVA EKIQR

Expression System E. coli

Purity > 85% by SDS-PAGE.


Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX67950-pro Image

3μg Human ZFAND5 protein (GTX67950-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).