

Human PAX8 protein, His tag

Cat. No. GTX67951-pro

Species	Human
---------	-------

Package

50 µg

Applications

Observed MW (kDa)	48 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	PBS, 20% Glycerol, 1mM DTT
--------	----------------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	---

Region/Sequence	Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPHNSIR SGHGGLNQLG GAFVNGRPLP EVVRQRIVDL AHQGVPRCDI SRQLRVSHGC VSKILGRYYE TGSIRPGVIG GSKPKVATPK VVEKIGDYKR QNPTMFAWEI RDRLLAEGVC DNDTVPSVSS INRIIRTKVQ QPFNLPMDS C VATKSLSPGH TLIPSSAVTP PES PQSDSLG STYSINGLLG IAQPGSDKRK MDDSDQDSCR LSIDSQSSSS GPRKHLRTDA FSQHHLEPLE C PFERQHYPE AYASPSHTKG EQEVNTLAMP MATPPTPPTA RPGASPTAC
-----------------	---

Expression System	E. coli
-------------------	---------

Purity	> 90% by SDS-PAGE.
--------	--------------------


Conjugation	Unconjugated
-------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67951-pro Image

3μg Human PAX8 protein (GTX67951-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).