

Human NELF-E protein, His tag

Cat. No. GTX67957-pro

Species Human**Package**

20 µg

Applications

Observed MW (kDa) 43 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 5mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMLVIPPG LSEEEALQK KFNKLKKKKK
ALLALKKQSS SSTTSQGGVK RSLSEQPVMD TATATEQAKQ LVKSGAISAI KAETKNSGFK RSRTLEGKLG DPEKGPVPTF
QPFQRSISAD DDLQESSRRP QRKSLYESFV SSSDRLRELG PDGEEAEGPG AGDGPPRSFD WGYEERSGAH SSASPPRSRS
RDRSHERNRD RDRDRERDRD RDRDRDRERD RDRDRDRDRD RERDRDRERD RDRDREGPFR RSDSFPERRA PRKGNTLYVY
GEDMTPTLLR GAFSPFGNII DLSMDPPRNC AFVTYKMEs ADQAVAEING TQVESVQLKV NIARKQPMLE AATGKSVWGS
LAVQNSPKGC HRDKRTQIVY SDDVYKENLV DGF

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67957-pro Image

3μg Human NELF-E protein (GTX67957-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).