

Human PACT protein, His tag

Cat. No. GTX67982-pro

Species Human

Package

50 µg

Applications

Observed MW (kDa) 34 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 2mM DTT, 1mM EDTA

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMQSRHR AEAPPLERED SGTFSLGKMI TAKPGKTPIQ VLHEYGMKTK NIPVYECERS DVQIHVPTFT FRVTVDITC TEGTSSKKLA KHRAAEAAIN ILKANASICF AVPDPLMPDP SKQPKNQ LNP IGSLQELAIH HGWRLPEYTL SQEGGPAHKR EYTTICRLES FMETGKGASK KQAKRNAAEK FLAKFSNISP ENHISLTNVV GHSLGCTWHS LRNSPGEKIN LLKRSLLSIP NTDYIQLLSE IAKEQGFNIT YLDIDELSAN GQYQCLAELS TSPITVCHGS GISCNAQSD AAHNALQYLK IIAERK

Expression System E. coli

Purity > 80% by SDS-PAGE.


Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX67982-pro Image

3μg Human PACT protein (GTX67982-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).