

## Human RFXANK protein, His tag

Cat. No. GTX67986-pro

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

Package

20 µg

## Applications

|                   |         |
|-------------------|---------|
| Observed MW (kDa) | 28 kDa. |
|-------------------|---------|

## Properties

|      |        |
|------|--------|
| Form | Liquid |
|------|--------|

|        |                                         |
|--------|-----------------------------------------|
| Buffer | 20mM Tris-HCl, 150mM NaCl, 20% Glycerol |
|--------|-----------------------------------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                                                                                            |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                                                                             |
|---------------|-----------------------------------------------------------------------------|
| Concentration | 0.25 mg/ml (Please refer to the vial label for the specific concentration.) |
|---------------|-----------------------------------------------------------------------------|

## Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMELTQPA EDLIQTQQTP ASELGDPEDP  
GEEAADGSDT VVLSLPCTP EPVNPEPDAS VSSPQGSSLK HSTTLNTRQR GNEVSALPAT LDCDNLVNKP DERGFTPLIW  
ASAFGEIETV RFLLEWGADP HILAKERESA LSLASTGGYT DIVGLLLERD VDINIYDWNG GTPLLYAVRG NHVKCVEALL ARGADLTTEA  
DSGYTPMDLA VALGYRKVQQ VIENHILKLF QSNLVPADPE

|                   |         |
|-------------------|---------|
| Expression System | E. coli |
|-------------------|---------|

|        |                    |
|--------|--------------------|
| Purity | > 85% by SDS-PAGE. |
|--------|--------------------|

|             |              |
|-------------|--------------|
| Conjugation | Unconjugated |
|-------------|--------------|

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX67986-pro Image

3μg Human RFXANK protein (GTX67986-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).