

Human APPBP1 protein, His tag

Cat. No. GTX68013-pro

Species	Human
---------	-------

Package

50 µg

Applications

Observed MW (kDa)	60 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 200mM NaCl, 10% Glycerol, 2mM DTT
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMALGKL LKEQKYDRQL RLWGDHGQEA
LESAHVCLIN ATATGTEILK NLVLPGIGSF TIIDGNQVSG EDAGNNFFLQ RSSIGKNRAE AAMEFLQELN SDVSGSFVEE
SPENLLDNDP SFFCRFTVVV ATQLPESTSL RLADVLWNSQ IPLLICRTYG LVGYMRIIK EHPVIESHPD NALEDLRLDK PFPELREHFQ
SYDLDHMEKK DHSHTPWIVI IAKYLAQWYS ETNGRIPTY KEKEDFRDLI RQGILKNENG APEDEENFEE AIKNVNTALN TTQIPSSIED
IFNDDRCINI TKQTPSFIL ARALKEFVAK EGQGNLPVRG TIPDMIADSG KYIKLQNVYR EKAKKDAAV GNHVAKLLQS
IGQAPESISE KELKLLCSNS AFLRVRCRS LAEEYGLDTI NKDEIISMD NPDNEIVLYL MLRAVDRFHK QQGRYPGVSN YQVEEDIGKL
KSCLTGLQE YGLSVMVKDD YVHEFCRYGA AEPHTIAAFL GGAAAEQVIK IITKQVIFN NTYIYSGMSQ TSATFQL

Expression System	E. coli
-------------------	---------

Purity	> 90% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68013-pro Image

3μg Human APPBP1 protein (GTX68013-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).