

Human PSTPIP1 protein, His tag

Cat. No. GTX68024-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 48 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPQLQF KDAFWCRDFT AHTGYEVLLQ
RLDGRKMCK DMEELLRQRA QAEERYGKEL VQIARKAGGQ TEINSLRASL DSLKQQMENV GSSHIQLALT LREELRSLEE
FRERQKEQRK KYEAVMDRVQ KSKLSLYKKA MESKKTIEQK CRDADDAEQA FERISANGHQ KQVEKSQNK A RQCKDSATEA
ERVYRQSI AQ LEKVRAWEQ EHRITCEAFQ LQEFDRILT RNALWVHSNQ LSMQCVKDDE LYEEVRLTLE GCSIDADIDS
FIQAKSTGTE PPAPVPYQNY YDREVTPLTS SPGIQPCSGM IKRFSGLLHG SPKTTSLAAS AASTETLTPT PERNEGVT A IAVQEIQGNP
ASPAQEYRAL YDYTAQNPDE LDLSAGDILE VILEGEDGWW TVERNGQRGF VPGSYLEKL

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68024-pro Image

3μg Human PSTPIP1 protein (GTX68024-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).