

Human Cyclin B2 protein, His tag

Cat. No. GTX68030-pro

Species Human

Package
50 µg

Applications

Observed MW (kDa) 45 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 5mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSHMALLRR PTVSSDLENI DTGVNSKVKS
HVTIRRTVLE EIGNRVTTTRA AQVAKKAQNT KVPVQPTKTT NVNKQLKPTA SVKPVQMEKL APKGPSPTPE DVSMKEENLC
QAFSDALLCK IEDIDNEDWE NPQLCSDYVK DIYQYLRQLE VLQSINPHFL DGRDINGRMR AILVDWLQVQV HSKFRLLQET
LYMCVGIMDR FLQVQPVSRL KLQLVGITAL LLASKYEEMF SPNIEDFVYI TDNAYTSSQI REMETLILKE LKFELGRPLP LHFLRRASKA
GEVDVEQHTL AKYLMELTLI DYDMVHYHPS KVAASASCLS QKVLGQGGKWN LKQQYYTGYT ENEVLEVMQH MAKNVVKVNE
NLTKFIAIKN KYASSKLLKI SMIPQLNSKA VKDLASPLIG RS

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX68030-pro Image

3μg Human Cyclin B2 protein (GTX68030-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).