

Human GRP1 protein, His tag

Cat. No. GTX68039-pro

Species Human**Package**

20 µg

Applications

Observed MW (kDa) 46 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMDEDGGG EGGGVPELDS LEEREELLDI
RRRKELIDD IERLKYEIAE VMTEIDNLTS VEESKTTQRN KQIAMGRKKF NMDPKKGIQF LIENDLLQSS PEDVAQFLYK GEGLNKTVIG
DYLGERDEFN IKVLQAFVEL HEFADLNLVQ ALRQFLWSFR LPGEAQKIDR MMEAFASRYC LCNPGVFQST DTCYVLSFAI
IMLNTSLHNH NVRDKPTAER FIAMNRGINE GGDLPCELLR NLYESIKNEP FKIPEDDGND LHTTFFNPDR EGWLLKLGGGR
VKTWKRWWFI LTDNCLYYFE YTTDKPRGI IPLENLSIRE VEDPRKPNCF ELYNPSHKGQ VIKACKTEAD GRVVEGNHVV YRISAPSPEE
KEEWMKSIKA SISRDPFYDM LATRKRRRIAN KK

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68039-pro Image

3μg Human GRP1 protein (GTX68039-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).