

Human Homer1 protein, His tag

Cat. No. GTX68063-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 40 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMGE QPIFSTRAHV
FQIDPNTKKK WVPTSKHAVT VSYFYDSTRN VYRIISLDGS KAIINSTITP NMTFTKTSQK FGQWADSRAN TVYGLGFSSE
HHLSKFAEKF QEFKEAARLA KEKSQEKMEI TSTPSQESAG GDLQSLTPE SINGTDDERT PDVTQNSEPR AEPTQNALPF
SHSSAISKHW EAEATLKGK NAKLTAALLE STANVKQWKQ QLAAYQEEAE RLHKRVTELE CVSSQANAVH THKTELNQTI
QELETCLKK EEEIERLKQE IDNARELQEQ RDSLTQKLQE VEIRNKDLEG QLSDLERLE KSQNEQEAFL NNLKTLLEIL DGKIFELTEL
RDNLAKLLEC S

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX68063-pro Image

3μg Human Homer1 protein (GTX68063-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).