

Human SCRN1 protein, His tag

Cat. No. GTX68092-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 46 kDa.

Properties

Form Liquid

Buffer PBS, 10% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMAAAPP YCFVAFPPRA KDGLVFGKN SARPRDEVQE VYFSAADHE PESKVECTYI SIDQVPRTYA IMISRPALW GAEMGANEHG VCIANEAIN REPAAEIAL LGMDLVRLGL ERGETAKEAL DVIVSLEEH GQGNYFEDA NSCHSFQSAY LIVDRDEAWV LETIGKYWAA EKVTEGVRCI CSQSLTTKM DAEHPELSY AQSQGWWTGE GEFNFSEVFS PVEDHLCGA GKDSLEKQEE SITVQTMNT LRDKASGVCI DSEFFLTAS GVSVPQNRSPCIHYFTGT PDPSRSIFKP FIFVDDVKLV PKTQSPCFGD DDPAPKEPRF QEKPDRRHEL YKAHEWARAI IESDQEQGRK LRSTMLELEK QGLEAMEEIL TSSEPLDPAE VGDLYDCVD TEIKFFK

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68092-pro Image

3μg Human SCR1 protein (GTX68092-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).