

Human G3BP2 protein, His tag

Cat. No. GTX68093-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 54 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 200mM NaCl, 30% Glycerol, 2mM DTT, 0.1mM PMSF

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSHMVMKEP SPLLVGREFV RQYYTLNKA
PEYLHRYFYGR NSSYVHGGVD ASGKPQEAHV GQNDIHHKVL SLNFSECHTK IRHVDAHATL SDGVVQVMG LLSNSGQPER
KFMQTFVLAP EGSVPNKFYV HNDMFRYEDE VFGDSEPELD ESEDEVEEEE QEERQPSPEP VQENANSQGY EAHPVTNGIE
EPLEESSHEP EPEPESETKT EELKPQVEEK NLEELKEEST TPPPAEPVSL PQEPPKPRVE AKPEVQSQPP RVREQRPRER PGFPPRGPRP
GRGDMEQNDS DNRRIIRYPD SHQLFVGNLP HDIDENELKE FFMSFGNVVE LRINTKGVGG KLPNFGFVVF DDSEPVRIL
IAKPIMFRGE VRLNVEEKKT RAARERETR GGD DRRDIRR NDRGPGGPRG IVGGGMMRDR DGRGPPPRGG MAQKLGSGRG
TGQMEGRFTG QRR

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68093-pro Image

3μg Human G3BP2 protein (GTX68093-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).