

Human FARSB protein, His tag

Cat. No. GTX68109-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 66 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPVSVK RDLLFQALGR TYTDEEFDEL
CFEFGLELDE ITSEKEIISK EQGNVKAAGA SDVVLKYIDV PANRYDLLCL EGLVRGLQVF KERIKAPVYK RVMPDGKIQK LIITEETAKI
RPFAVAVALR NIKFTKDRYD SFIELQEKHL QNICRKRALV AIGTHDLDTL SGPFTYTAKR PSDIKFKPLN KTKEYTACEL MNIYKTDNHL
KHYLHIIENK PLYPVIYDSN GVLSMPPII NGDHSRITVN TRNIFIETG TDFTKAKIVL DIIVTMFSEY CENQFTVEAA EVVFPNGKSH
TFPELAYRKE MVRADLINKK VGIRETPENL AKLLTRMYLK SEVIGDGNQI EIEIPPTRAD IIAACDIVED AAIAYGYNNI QMTLPKTYTI
ANQFPLNKL ELLRHDMAAA GFTEALTFAL CSQEDIADKL GVDISATKAV HISNPKTAEF QVARTLLPG LLKTIAANRK MPLPLKLEI
SDIVIKDSNT DVGAKNYRHL CAVYYNKNPG FEIIGHLLDR IMQLLDVPPG EDKGGYVIKA SEGPAFFPGR CAEIFARGQS
VGKLGVLHPD VITKFELTMP CSSLEINVGP FL

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68109-pro Image

3μg Human FARSF protein (GTX68109-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).