

Human G3BP1 protein, His tag

Cat. No. GTX68122-pro

Species	Human
---------	-------

References (2)

Package

50 µg

Applications

Observed MW (kDa)	52 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 2mM DTT
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, C-terminal His-Tag: MVMEKPSPLL VGREFVRQYY TLLNQAPDML HRFYGNSSY VHGGGLDSNGK
PADAVYGQKE IHRKVMSQNF TNCHTKIRHV DAHATLNDGV VQVMGLLSN NNQALRRFMQ TFVLAPEGSV ANKFYVHNDI
FRYQDEVFGG FVTEPQEES EEEVEEERQ QTPEVVPDDS GTFYDQAVVS NDMEEHLEEP VAEPEPDPEP EPEQEPVSEI QEEKPEPVLE
ETAPEDAQKS SSPAPADIAQ TVQEDLRTFS WASVTSKNLP PSGAVPVTGI PPHVVKVPAS QRPESKPES QIPPQRPQRD
QRVREQRINI PPQRGPRPIR EAGEQGDIEP RRMVVRHPDSH QLFIGNLPHE VDKSELKDFE QSYGNVVELR INSGGKLPNF
GFVVFDDSEP VQKVLNRPI MFRGEVRLNV EEKKTRAARE GDRRDNRRLRG PGGPRGGLGG GMRGPPRGGM VQKPGFGVGR
GLAPRQVEHH HHHH

Expression System	E. coli
-------------------	---------

Purity	> 85% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68122-pro Image

3μg Human G3BP1 protein (GTX68122-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).