

Human GLYAT protein, His tag

Cat. No. GTX68132-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 34 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MMLPLQGAQM LQMLEKSLRK SLPASLKVYG TVFHINHGPN FNLKAVVDKW PDFNTVVVCP QEQDMTDDL HYNTYQIYS KDPQNCQFL GSPGLINWKQ HLQIQSSQPS LNEAIQNLA IKSFKVKQTQ RILYMAAETA KELTPFLKS KILSPSGGKP KAINQEMFKL SSMDVTHAHL VNKFVHFHGGN ERSQRFIERC IQTFPTCLL GPEGTPVCWD LMDQTGEMRM AGTLPEYRLH GLVTYVIYSH AQLGKLGFP VYSHVDYSNE AMQKMSYTLQ HVPIPRSWNQ WNCVPL

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68132-pro Image

3μg Human GLYAT protein (GTX68132-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).