

Human BPNT1 protein, His tag

Cat. No. GTX68146-pro
Species Human

Package

50 µg

Applications
Observed MW (kDa) 33 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 5mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMASST NTVMRLVAS
 AYSIAQKAGM IVRRVIAEGD LGIVEKTCAT DLQTKADRLA QMSICSSLAR KFPKLTIGE EDLPSEEDVQ ELIEDSQWEE ILKQPCPSQY
 SAIKEEDLVV WVDPLDGTKE YTEGLLDNVT VLGIAIYEGK AIAGVINQPY YNYEAGPDAV LGRTIWGVLG LGAFGFQLKE VPAGKHIITT
 TRSHSNKLVV DCVAAMNPDA VLRVGGAGNK IIQLIEGKAS AYVFASPGCK KWDTCPEVI LHAVGGKLTG IHGNNVLQYHK
 DVKHMNSAGV LATLRNYDYY ASRVPEIKN ALVP

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68146-pro Image

3μg Human BPNT1 protein (GTX68146-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).