

Human beta Tubulin 3/ Tuj1 protein, His tag

Cat. No. GTX68147-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 50 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMREIVHI QAGQCGNQIG AKFWEVISDE
HGIDPSGNYV GDSDLQLERI SVYYNEASSH KYVPRAILVD LEPGTMDSVR SGAFGHLFRP DNFIGQSGA GNNWAKGHYT
EGAELVDSVL DWRKECENC DCLQGFQLTH SLGGGTGSGM GTLLISKVRE EYPDRIMNTF SWPSPKVSD TVPEPYNATL
SIHQLVENTD ETYCIDNEAL YDICFRTLKL ATPYGDLNH LVSATMSGVT TSLRFPQQLN ADLRKLAVNM VPFRLHFFM
PGFAPLTARG SQQYRALVTP ELTQQMFDAK NMMAACDPRH GRYLTVATVF RGRMSMKEVD EQMLAIQSKN SSYFVEWIPN
NVKVAVCDIP PRGLKMSSTF IGNSTAIQEL FKRISEQFTA MFRKAFLEHW YTGEGMDEME FTEAESNMND LVSEYQQYQD
ATAEEEGEMY EDDEEESEAQ GPK

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68147-pro Image

3μg Human beta Tubulin 3/ Tuj1 protein (GTX68147-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).