

Human HAX1 protein, His tag

Cat. No. GTX68154-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

Applications

Observed MW (kDa)	32 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 100mM NaCl, 10% Glycerol
--------	---

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MSLFDLFRGF FGFGPGRSHR DPFFGGMTRD EDDDEEEEE GGSWGRGNPR FHSPQHPPPEE FGFGFSFSPG GGIRFHDNFG FDDLVRDFNS IFSDMGAWTL PSHPPPELPGP ESETPGERLR EGQTLRDSML KYPDSHQPRI FGGVLES DAR SESPQAPDW GSQRPFHREF DVWPMDPHPR TREDNDLDSQ VSQEG LGPVL QPQPKSYFKS ISVTKITKPD GIVEERRTV DSEGR TETTV TRHEADSSPR GDPESRPPA LDDAFSILD L FLGRWFRSR
-----------------	---

Expression System	E. coli
-------------------	---------

Purity	> 80% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------


Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX68154-pro Image

3μg Human HAX1 protein (GTX68154-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).